



## **MiniSim 1000**

**Multi Parameter Patient Simulator**



## **INSTRUCTION MANUAL**

110 Toledo Street  
Farmingdale N.Y. 11735  
Phone: 631-531-0100 Fax: 631-531-0101  
[www.Netech.org](http://www.Netech.org)

## Copyright

Copyright © 2003 by Netech Corporation. All rights reserved. No part of this publication may be reproduced or transmitted in any form other than for the purchaser's personal use without written permission from Netech Corporation.

## Quality Assurance

Netech Corporation is ISO 9001-2000 Registered. This instrument was thoroughly tested and inspected according to Netech's ISO 9001-2000 quality standards (ISO/IEC 17025) and test procedures, and was found to meet those specifications when it was shipped from the factory.

## Warranty

Netech warrants the MiniSim 1000 Multi-Parameter Patient Simulator against defects in materials and workmanship for one year from the date of original purchase.

During the warranty period, we will repair or, at our option, replace at no charge a product that proves to be defective, provided you return the product shipping prepaid to Netech Corporation. This warranty does not apply if the product has been damaged by accident or misuse, or as the result of service or modification by other than Netech Corporation, or if its serial number is defaced or removed.

Netech reserves the right to discontinue the MiniSim 1000 at any time, and change its specifications, price, or design without notice and without incurring any obligation. Netech guarantees availability of service parts for 5 years after the manufacture of the unit is discontinued.

The warranty is void if you elect to have the unit serviced and / or calibrated by someone other than Netech.

***The purchaser assumes all liability for any damages or bodily injury that may result from the use or misuse of the unit by the purchaser, his employees, agents or customers.***

**In no event shall Netech Corporation be liable for consequential damages.**

## **Trademarks**

Netech and MiniSim 1000 are trademarks of Netech Corporation. Any other trademark names used in this manual are only for editorial purposes and the benefit of the respective trademark owner, with no intention of improperly using that trademark.

## Table of Contents

---

### 1. General Overview

Introduction	-----	06
Specifications	-----	08
Accessories	-----	10

### 2. Operating Instructions

Controls and Indicators	-----	08
Baseline ECG Waveform	-----	09
Performance Waveforms	-----	11
Arrhythmia Waveforms	-----	12
Automatic Testing Sequences	-----	13
Respiration Waveforms/ Apnea	-----	16
Blood Pressure waveforms	-----	16
Temperature	-----	17
Calibration	-----	17

## GENERAL OVERVIEW

---

### 1.1 Introduction

#### GENERAL DESCRIPTION

The MiniSim 1000 Multi-Parameter Patient Simulator is an advanced microcontroller based instrument. It is designed to simulate patient signals of ECG, arrhythmia, invasive blood pressure, respiration, and temperature. The device also simulates square, sine, triangle, and pulse performance waveforms.

The easy to operate MiniSim 1000 is menu driven via eight tactile feel soft keys. All functions are displayed on a two line sixteen character LCD display.

The small hand held instrument is powered by one 9 Volt battery or an optional AC adapter. The MiniSim 1000 is CE marked and shipped with a Certificate of Calibration traceable to the NIST.

The MiniSim 1000 is a rugged instrument that performs its simulations quickly, accurately, and with ease.

### 1.2 Specifications

**ECG:** 12 leads with independent outputs referenced to RL.

#### NORMAL SINUS RHYTHM (NSR)

**ECG Rates:** 30, 60, 70, 80, 90, 100, 120, 150, 180, 210, 240, 270, 300, and 350 BPM.  
Accuracy 0.5%.

**Amplitudes:** 0.15, 0.3, 0.5, 1.0, 2.0, 3.0, 4.0, and 5.0 mV on Lead II. Lead I is 0.6 X Lead II, Lead III is 0.4 X Lead II, and Lead V is the same as Lead II.

**High Level:** 500 X low level output on Lead II.  
Accuracy: 2% (1-5mV).

#### PERFORMANCE WAVEFORMS (PERF)

##### Sine, Square, Triangle, Pulse

**Frequencies:** 0.1 to 0.9 in 0.1 Hz steps.  
1.0 to 9.0 in 1.0 Hz steps.  
10 to 100 in 10 Hz steps.

**Accuracy:** 1%

**Amplitudes:** 0.1, 0.2, 0.5, 1.0, 2.0, 3.0, 4.0, and 5.0 mV on Lead II.  
**Pulse:** 20 ms pulse of 1mV amplitude repeated at 4 second intervals.

#### R WAVE DETECTION (RWD)

**Widths:** 10, 40, 50, 60, 70, 80, 90, 100, 110, and 120 ms.  
**Amplitudes:** 0.15, 0.3, 0.5, 1.0, 2.0, 3.0, 4.0, 5.0 mV on Lead II.

### **PACER (PCR)**

**AP:** Atrial Pacer set at 70 BPM.  
**ASP:** Asynchronous Pacer – Ineffective Pacing.  
**Pacer Amplitudes:** + 2 mV and –2mV.  
**QRS Amplitudes:** 1 mV  
**PCR:** Pacer pulses alone.  
**Pacer Widths:** 0.1, 0.2, 0.5, 1.0, 1.5, and 2.0 ms.  
**Amplitudes:** 0.15, 0.3, 0.5, 1.0, 2.0, 3.0, 4.0, and 5.0 mV.  
**VP:** Ventricular Pacer set at 70 BPM.  
**AVP:** Atrial Ventricular Pacer set at 70 BPM.  
**QRS Amplitude:** 1 mV.  
**Pacer Amplitude:** -2mV.

### **RESPIRATION**

**Baseline Impedance:** 250, 500, 750, and 1000 Ohms.  
**Delta Impedance:** 0.1, 0.5, 1.0, and 1.5 Ohms.  
**Accuracy:** 10%.  
**Rates:** 15, 30, 60, and 120 BPM.  
**Accuracy:** 1%  
**Apnea:** Off, continuous, 12 seconds and 32 seconds.  
**Lead Configuration:** Leads I and II.

### **BLOOD PRESSURE**

**Impedance:** 350 Ohms.  
**Excitation:** 2 to 16 Volts.  
**Sensitivity:** 5 $\mu$ V/V/mmHg.  
**Static:** 0, 5, 10, 20, 25, 30, 40, 50, 100, 150, 200, and 300 mmHg.  
**Dynamic:** 50/10, 60/20, 70/30, 80/40, 100/60, and 120/80 mmHg. Dynamic waveforms track all NSR rates.  
**Accuracy:** 1%

### **TEMPERATURE**

**Compatibility:** YSI 400 and 700 Series.  
**Temperature:** 25 and 37 degrees Centigrade.  
**Accuracy:** 2% of setting

### **POWER REQUIRMENTS:**

One 9 Volt alkaline battery or optional AC adapter.

### **PHYSICAL CHARACTERISTICS:**

**Dimensions:** 5.5 X 3.5 X 1.5 inches (13.9 X 8.9 X 3.8 cm).  
**Weight:** 10 oz ( 0.3 kg).

### **TEMPERATURE REQUIREMENTS:**

**Operating:** 59 to 95° F (15 to 35°C).  
**Storage:** 32 to 131° F (0 to 55°C).

### **CALIBRATED DC OUTPUTS:**

10 Calibrated DC outputs from 0.04 to 2 mV.

## OUTPUT CONNECTIONS:

**Part # 303**, Open Ended Pressure cable:

Pressure: Single Pressure  
Pin 1 + Excitation  
Pin 4 + Signal  
Pin 2 - Signal  
Pin 5 - Excitation

Pressure: Dual Pressure  
Pin 3 + Signal  
Pin 2 - Signal

**Part # 304**, Open Ended Temperature Cable:  
3.5 mm Stereo Jack

YSI 400 Tip, Ring, Shield.  
YSI 700 Tip, Ring.

High Level Output: 3.5 mm Stereo Jack

Signal Positive Tip  
Signal Ground Shield  
Input/UnusedRing

**Note:** Specifications are subject to change without notice.

### 1.3 ACCESSORIES

Description	Part Number
User Manual	300-Manual
Pressure Cable (Open Ended)	303
Temperature Cable (Open Ended)	304

### OPTIONAL ACCESSORIES

Description	Part Number
Soft Carrying Case	301
110 VAC Adapter	302
220 VAC Adapter	302-220
ECG Snap Studs	1000
Space Labs Dual	305
Space Labs Single	305-S
Datascope Dual	306
Datascope Single	310
Hewlett Packard Dual	307
Hewlett Packard Single	308
Fukuda Denshi Dual	309
Bard Dual	311
Bard Single	311-S
Mennen Medical Single	312
Ivy Biomedical Single	314

**Other Interface cables available. Please call or send an email for more information.**

**(800) 547-6557**  [Sales@Netech.org](mailto:Sales@Netech.org)

## 2. OPERATING INSTRUCTIONS

### Menu Descriptions

The following are the menu sequence displays, their abbreviations, definitions, and default settings.

SIMULATOR MENU			
Ecg	Resp	Bp	Temp

**Ecg:** ECG

**Default Setting:** Normal Sinus Rhythm (NSR) at 80 BPM with Amplitude of 1 mV

**Resp:** Respiration

**Default Setting:** 30 BPM, Impedance 500 Ohms, Delta R 1.0 Ohm, and Apnea Off

**Bp:** Blood Pressure


**Default Setting:** 120/80 Dynamic Setting  
0 Static Setting

**Temp:** Temperature

**Default Setting:** YSI 700 series set to 25°C  
YSI 400 series set to 25°C

## CONTROLS AND INDICATORS

1. On-Off Key
2. ECG Lead Snaps
3. LCD Display
4. Temperature Connector
5. Blood Pressure Connector
6. Hi-Level ECG Output
7. Function Keys: 'F1', 'F2', 'F3', 'F4'.
8. Arrow Keys: '<' and '>'.
9. Menu Key



The operating menu of the MiniSim 1000 is arranged in a tree structure. When turned on, the microcontroller initializes a test routine and displays the model and software version numbers and then switches to the Main Menu.

During operation, the 'MENU' key returns the currently displayed menu to the previous menu.



The 'F1', 'F2', 'F3', and 'F4' keys select the parameter for simulation, the specific type of parameter, and the functional characteristics of the parameter.

The '<' and '>' arrow keys move through the available choices under the selected parameter.

**Basic Operating Instructions:**

1. Connect the MiniSim 1000 to a patient monitor using the ECG lead snaps, correct blood pressure interface cable, and correct temperature interface cable.
2. Turn the Minisim 1000 On.
3. Select the parameter for simulation, 'ECG', 'Resp', 'BP', or 'Temp'.
4. Continue to select choices under the parameter to be simulated until the functional characteristics have been chosen.
5. Each Parameter has default settings.

<b>ECG Menu</b>			
<b>Base</b>	<b>Perf</b>	<b>Arth</b>	<b>Aut</b>

**Base:** Baseline ECG  
**Perf:** Performance Waveforms  
**Arth:** Arrhythmia Waveforms  
**Aut:** Automatic Test Sequences

<b>Baseline ECG</b>			
<b>NSR</b>	<b>PCR</b>	<b>ST</b>	<b>RWD</b>

**NSR:** Normal Sinus Rhythm  
**PCR:** Pacer Waveforms  
**ST:** ST Segment Analysis Waveforms  
**RWD:** R Wave Detection

<b>NSR/</b>	
<b>Rate</b>	<b>Ampl</b>

Rate and Amplitude selections under NSR are:

**Rate:** 30 60 70 80 90 100 120 150 180 210 240 270 300 350  
**Ampl:** .15 0.3 0.5 1.0 2.0 3.0 4.0 5.0

<b>Pacer Waves</b>	
<b>Atr</b>	<b>Vent</b>

Selecting **Atr** branches to:

**Atr\_Pcr Waves**  
**+ve AP PCR ASP**

**+ve:** Changes pacer pulse from positive (+ve) to negative (–ve).

**AP:** Atrial Pacer – Normal Paced Rhythm

QRS with rate of 60 BPM and Pacemaker pulses with amplitudes of + 2 mV and –2 mV with duration of 0.1 ms or 2 ms with a normally paced QRS T (QRS amplitude of 1 mV, duration of 100 ms, T wave amplitude of 0.2 mV, duration of 180 ms, and Q-T interval of 350ms)

**PCR:** Pacer Pulses Alone. Amplitude default setting of 1.0 mV and width of 1.0 ms.

**Widt:** width selections

0.1 0.2 0.5 1.0 1.5 2.0

**Ampl:** amplitude selections

.15 0.3 0.5 1.0 2.0 3.0 4.0 5.0

**ASP:** Asynchronous Pacer – Ineffective pacing.

A non-synchronized waveform that combines QRS waves at 30 BPM and other specifications as in AP above with pacer waves with a 0.1 ms or 2 ms width and amplitude of  $\pm 2$  mV at 80 BPM.

Note: 0.1 ms pacer width for AP and ASP can only be selected by going to 'PCR'. A selection other than 0.1 ms produces 2 ms wide pulses.

Selecting **Vent** branches to:

**Vtr\_Pcr Waves**  
**VP AVP**

**VP:** Ventricular Pacer at 70 BPM

**AVP:** Atrial Ventricular Pacer at 70 BPM

Selecting **ST** branches to:

**ST-SEGMENT**  
**ele dep MI TalT**

ST segment analysis waveforms are divided into four classes:

**ele:** ST Elevation

**dep:** ST Depression

**MI:** Myocardial Infarction

**TalT:** Tall T Wave Rejection

**ST\_ELAVATION**  
**7% 13% 20% Flat**

**ST Elevation:** 7, 13, and 20% DC levels of ST Elevation may be selected. Example: at a QRS amplitude of 1 mV, ST segments are produced at positive DC levels of 70, 130, and 200 micro volts.

Waveforms may be selected with a 'Flat', a positive '+sl', or a negative '-sl' slope. The ST segment is proportional to the ECG amplitude setting.

**ST\_DPRESION**  
**7% 13% 20% Flat**

**ST Depression:** Identical to ST Elevation except the ST segments are now depressed.

Selecting **MI Myocardial Infarction** branches to:

<b>Myocard_Inf</b>			
<b>Isc</b>	<b>Inj</b>	<b>Inf</b>	<b>linf</b>

**Isc: Ischemia**  
**Inj: Injury**  
**Inf: Infarction**  
**linf: Inferior Infarction**

**Isc:** Ischemia : Normal Sinus Rhythm (NSR) with fully inverted T waves. This is a condition of reduced blood supply to the heart in a normal patient.

**INJ:** Injury. A waveform with ST elevation of 20% with a negative slope and inverted T wave.

**Inf:** Infarction. Normal Sinus Rhythm (NSR) with a large Q wave with the amplitude increased six times and the width increased three times compared to normal.

**linf:** Inferior Infarction. A waveform with the Q wave modified as in Infarction and the ST segment elevated 7% as in ST Elevation.

Selecting **TaIT** branches to:

<b>TaIT/Rate = 80</b>	
<b>Rate</b>	<b>Ampl</b>

An 80 BPM QRS test signal of 1 mV amplitude and 100 ms duration is generated with a T wave duration of 180 ms and Q-T interval of 350 ms.

The T wave amplitude may be varied from 0 to 1.2 mV in steps of 0.1 mV .

**Rate:** 80 BPM  
**Amplitude:** 0.0 0.1 0.2 0.3 0.5 0.6 0.7 0.8 0.9 1.0 1.1 1.2

Selecting **RWD** branches to:

<b>RWD/</b>	
<b>Widt</b>	<b>Ampl</b>

**RWD:** R Wave Detection. A QRS waveform is generated at 70 BPM with selectable width and amplitude changes.

**Widt:** Width default setting 100 msec  
 10 40 50 60 70 80 90 100 110 120

**Ampl:** Amplitude default setting 1.0 mV  
 .15 0.3 0.5 1.0 2.0 3.0 4.0 5.0

<b>Perf. Waves</b>			
<b>SIN</b>	<b>SQR</b>	<b>TRI</b>	<b>PLS</b>

**SIN:** Sine Wave  
**SQR:** Square Wave  
**TRI:** Triangle Wave  
**PLS:** Pulse Wave

<b>SIN/</b>	
<b>Freq</b>	<b>Ampl</b>

<b>Freq</b>	<b>SQR/ Ampl</b>
<b>Freq</b>	<b>TRI/ Ampl</b>

The frequency and amplitude default settings and choices are the same for the sine, square and triangle performance waveforms.

**Freq:** Frequency default setting 1.0 Hz

0.9 1.0 2.0 3.0 4.0 5.0 6.0 7.0 8.0 9.0 10 20 30 40 50 60 70 80 90 100  
0.1 0.2 0.3 0.4 0.5 0.6 0.7 0.8

**Ampl:** Amplitude default setting 2.0 mV

2.0 3.0 4.0 5.0 0.2 0.5 1.0

**Pulse wave default setting:** A pulse wave is generated at 4 second intervals with an amplitude of 1mV and width of 20 ms.

Selecting **Arr** branches to:

<b>Arrhythmia Menu</b>		
<b>Atr</b>	<b>AC</b>	<b>Vent</b>

Arrhythmias are divided into Atrial, Atrial Conduction, and Ventricular waveforms. Normal Sinus Rhythm (NSR) at 80 BPM is the default waveform in this menu. Function keys select the desired arrhythmia and the MENU key clears the arrhythmia to NSR.

The following are the arrhythmia definitions:

#### **Atrial: Atrial Arrhythmias**

**SA: Sinus Arrhythmia:** The ECG rate uniformly increases and decreases continuously. The pattern is cyclic with rates changing in the following order: 60, 70, 80, 90, 100, 90, 80, 70, 60 BPM.

**M80: Missing Beat:** Normal Sinus Rhythm is generated at 80 BPM with every 10<sup>th</sup> beat missing.

**AFLT: Atrial Flutter:** Varying ECG rates with 12 cycles at 60 BPM for 12 seconds, 9 cycles at 90 BPM for 6 seconds, 15 cycles at 150 BPM for 6 seconds repeating with large P waves at 300 BPM. This corresponds to ventricular responses of 5:1 for 12 seconds, 3:1 for 6 seconds, and 2:1 for 6 seconds.

**AFB: Atrial Fibrillation:** Irregular QRS complexes with no P waves and constantly changing R-R intervals are generated. The rate varies in a cyclic fashion at 30, 60, 70, 80 and 30 BPM with low amplitude oscillations on the baseline.

**PAT: Paroxysmal Atrial Tachycardia:** NSR is generated at 180 BPM with inverted P waves.

**NODL: Junctional Premature Contraction:** NSR is generated at 80 BPM with a short PR interval. The QRS starts immediately following the P wave.

#### **AC: Atrial Conduction Arrhythmias**

**AB1: First Degree AV Block:** The QRS is generated at 80 BPM, the P wave precedes the QRS by a fixed but prolonged PR interval > 0.2 seconds ( PR interval = 0.26 seconds).

**MB1: Second Degree AV Block: Mobitz I: Wenckebach:** The QRS is generated at 80 BPM. There is a progressive lengthening of the PR interval with intermittent dropped beats. The PR intervals are 170, 230, and 310 ms.

**MB2: Second Degree AV Block: Mobitz II:** The QRS is generated at 80 BPM with every 4<sup>th</sup> QRS missing. The PR interval is constant at 170 ms.

**AB3: Third degree AV Block:** The P wave and QRS are independent of each other. The P wave is generated at 80 BPM and the QRS is generated at 50 BPM.

**RBB: Right Bundle Branch Block:** A prolonged QRS (>0.12 sec) is generated at 80 BPM. The resulting QRS looks like the letter "M".

**LBB: Left Bundle Branch Block:** A widened QRS is generated at 80 BPM with a large wide S wave.

**LAH: Left Anterior Hemiblock:** A QRS is generated at 80 BPM with an S wave larger than the R wave.

#### **Ventricular: Ventricular Arrhythmias**

**PV1: Premature Ventricular Contraction 1:** NSR is generated at 80 BPM. Each time the 'F1' function key is pressed one PVC is generated.

**PV3: Premature Ventricular Contraction 3:** NSR is generated at 80 BPM. Each time the 'F2' function key is pressed 3 PVCs are generated.

**PV6: Premature Ventricular Contraction 6:** NSR is generated at 80 BPM. Each time the 'F3' function key is pressed 6 PVCs are generated.

**PV12: Premature Ventricular Contraction 12:** NSR is generated at 80 BPM. Each time the 'F4' function key is pressed 12 PVCs are generated at different intervals.

**PV24: Premature Ventricular Contraction 24:** NSR is generated at 80 BPM. Each time the 'F1' function key is pressed 24 PVCs are generated at different intervals.

**BGY: Bigeminy:** NSR is generated at 80 BPM with every other beat as a PVC.

**TGY: Trigeminy:** NSR is generated at 80 BPM with every third beat as a PVC.

**PVC: Premature Ventricular Contraction:** Continuous PVCs are generated at 80 BPM.

**VFLT: Ventricular Flutter:** Sine waves at 240 BPM are generated with irregular amplitudes.

**VFB: Ventricular Fibrillation:** A totally irregular waveform is generated with chaotic undulations of the baseline.

**VTC: Ventricular Tachycardia:** A fast moving series of PVCs is generated at 210 BPM.

**PVR: Right Focal PVC:** NSR is generated at 80 BPM with every 10<sup>th</sup> beat a right focal PVC.

Selecting **Aut** branches to:

<b>&lt;Auto Sequence&gt;</b>			
<b>RWD</b>	<b>TalT</b>	<b>PPR</b>	<b>TAC</b>

The Automatic Test Sequence generates test waveforms in sequences according to AAMI requirements eliminating the need for the user to make numerous manual selections.

The automatic test sequences are:

**RWD: R Wave Detection**  
**TalT: Tall T Wave Rejection**  
**PPR: Pacemaker Pulse Rejection**  
**TAC: Time for Alarm for Tachycardia**

The automatic sequence is initiated when the test is selected and continues until all of the test patterns have been generated. During the test sequence the values generated are displayed on the LCD. At the completion of the automatic test, the display will return to its steady mode.

**RWD: R Wave Detection.** R waves are generated with three varying parameters of amplitude, width, and rate. The three parameter values are displayed on the LCD while generated. Each test waveform is displayed for 20 seconds.

The complete waveform test sequence is performed in three separate sets.

Set 1: Variable Parameters

Amplitude: 0.5, 2, 5 mV  
 Width: 100, 70, 120 ms  
 Rate: 80, 30, 210 BPM

Set 1 Complete Test Sequence

Amplitude (mV)	Width (ms)	Rate (BPM)	Time (sec)
0.5	100	80, 30, 210	20
	70	80, 30, 210	20
	120	80, 30, 210	20
2.0	100	80, 30, 210	20
	70	80, 30, 210	20
	120	80, 30, 210	20
5.0	100	80, 30, 210	20
	70	80, 30, 210	20
	120	80, 30, 210	20

As the test progresses, the indicated heart rate displayed on the patient monitor should be within  $\pm 10\%$  or  $\pm 5$  BPM whichever is greater of the applied rate.

Set 2: Variable Parameters

Amplitude: 0.15 mV  
 Width: 70, 120 ms  
 Rate: 30, 210 BPM

Set 2 Complete Test Sequence

Amplitude (mV)	Width (ms)	Rate (BPM)	Time (sec)
0.15	70	30, 210	20
0.15	120	30, 210	20

The monitor will not respond to the waveforms in this sequence Set.

Set 3: Variable Parameters

Amplitude: 1.0 mV  
 Width: 10 ms  
 Rate: 30, 210 BPM

Set 3 Complete Test Sequence

Amplitude (mV)	Width (ms)	Rate (BPM)	Time (sec)
1.0	10	30, 210	20

The patient monitor will not respond to the waveforms in the Set 3 test sequence.

**TalT: Tall T Wave Rejection.** QRS and T waves are generated with the following values:

QRS: Rate           80 BPM  
      Amplitude       1 mV  
      Width           100 ms  
      T Wave Duration 180 ms  
      QT Interval     350 ms  
T Wave Amplitude 0.0, 0.2, 0.4, 0.6, 0.8, 1.0, and 1.2 mV .

In the automatic test sequence the T Wave Amplitude steps through the changes at one minute intervals. The display indicates the T Wave Amplitude and the QRS Rate.

As the T wave amplitude increases, the first value at which the patient monitor counts the T wave at 80 ± 8 BPM should be noted. This value should match the patient monitor manufacturer's specification.

**PPR: Pacemaker Pulse Rejection.** The test sequence cycles through normal paced rhythm (AP), ineffective pacing (Asynchronous Pacing ASP), and pacemaker pulses alone (PCR).

For normal pacing (AP),QRS and pacer waves are generated with the following values:

QRS: Amplitude     1 mV  
      Width         100 ms  
T wave: Amplitude   0.2 mV  
      Duration      180 ms  
Q-T Interval        350 ms  
R-R Interval        1 Sec  
Pacer: Amplitude    2 mV, -2 mV  
      Width         2 ms, 0.1 ms

For ineffective pacing (ASP), the values of the QRS and pacer waves are the same as normal pacing except for the QRS rate that becomes 30 BPM and the pacer rate that becomes 80 BPM.

During the ASP and AP test sequences the display will show the pacer amplitude, the pacer width, and the QRS rate.

For pacemaker pulses alone (PCR) the values generated are:

Pacer: Rate         60 BPM  
      Width         2.0 ms, 0.1 ms  
      Amplitude     2 mV, -2 mV

During the PCR test sequence the display will show the pacer amplitude, the pacer width, and the QRS rate.

Each set of values in the test sequence is generated and displayed for 20 seconds.

**TAC: Time to Alarm for Tachycardia.** The TAC test is designed to measure the time it takes for the patient monitor to alarm after the onset of ventricular tachycardia. The low and high alarms on the patient monitor should be set at 60 BPM 100 BPM before starting the test.

This auto test sequence generates a QRS wave form at the rate of 80 BPM alternating with a ventricular tachycardia waveform with rates of 206 and 195 BPM and amplitudes of 1.0, 0.5, 2.0, and 4.0 mV.

Each waveform is generated for 20 seconds and the display will show the amplitude and the rate of the waveform being generated.

The following is the sequence of waveforms:

QRS		Ventricular Tachycardia		
Rate (BPM)	Amp (mV)	Width (ms)	Rate (BPM)	Amp (mV)
80	1	100	206	1.0
80	1	100	206	0.5
80	1	100	206	2.0

QRS		Ventricular Tachycardia		
Rate (BPM)	Amp (mV)	Width (ms)	Rate (BPM)	Amp (mV)
80	1	100	195	2.0
80	1	100	195	1.0
80	1	100	195	4.0

### RESPIRATION:

Select Resp in the Simulator Menu to access the Respiration selections.

Resp Menu			
Rate	Imp	dR	Apne

Respiration waveforms are generated with four selectable rates, baseline impedances, and delta impedance variations.

The Respiration default settings are:

Respiration Rate	30 BPM
Impedance	500 Ohms
Delta Impedance	1.0 Ohm
Apnea	Off

The value changes that may be made are:

Rate: 15, 30, 60, 120 BPM.

Impedance: 250, 500, 750, and 1000 Ohms.

Delta Impedance: 0.1, 0.5, 1.0, and 1.5 Ohms.

### Apnea:

Resp/Apne=Off			
Off	Cont	12s	32s

**Off:** Apnea is absent. Normal respiration waveforms are generated.

**Cont:** Continuous Apnea. No respiration waveforms are generated.

**12s:** No respiration waveform is generated for 12 seconds.

**32s:** No respiration waveform is generated for 32 seconds.

### BLOOD PRESSURE:

Select Bp in the Simulator Menu to access Blood Pressure selections.

BP1 Menu	
Dyna	Stat

Two blood pressure waveforms are generated with selections of 12 static and 6 dynamic values.

BP1 values are selected from the listed choices and BP2 values are ½ of those selected for BP1.

The Dynamic pressure waveforms track the Normal Sinus Rhythm rates.

The default settings for pressure values are:

BP1 Dynamic: 120/80 Static: 0

BP2 Dynamic: 60/40 Static: 0

The Dynamic pressure value selections are:

100/60      120/80

50/10      60/20

70/30      80/40

The Static pressure value selections are:

0   5   10   20


25   30   40   50

100   150   200   300



## TEMPERATURE:

Select Temp in the Simulator Menu to access Temperature value selections.



Temperature simulation is provided for both YSI 400 and YSI 700 standards. The temperature default setting for both is 25 degrees Centigrade.

The Temperature value selections are:

25 C 37 C for both YSI 400 and YSI 700.

## CALIBRATION

The MiniSim 1000 is shipped from the factory fully calibrated with a Certificate of Calibration traceable to the NIST. Annual calibration is recommended. Contact the Customer Service Department for instructions for returning the instrument to Netech for either calibration or repair.

### Extended Warranty Option:

If the instrument is returned to Netech Technical Service for its recommended annual calibration, the standard one year warranty is extended for a second year.

## Appendix A



ISO 9001 - 2000 Registered  
SBA 8 (a) Certified

**CATALOG 2007**

**Biomedical / Industrial**  
**Test & Measurement Instruments**

### Products

Arrhythmia Simulator  
Defibrillator Analyzer  
Pacemaker Analyzer  
Patient Simulator  
Safety Analyzer  
Pressure Meter  
Tachometer  
Flow Meter  
Service Kits  
SPO2 Simulator  
NIBP Simulator  
Ventilator Tester  
Laser Power Meter  
Ultrasound Watt Meter  
Light Meter



**Innovative Instruments Since 1987**



110 Toledo Street, Farmingdale, NY 11735  
Phone: 631-531-0100, Fax: 631-531-0101

Biomedical Site: <http://www.Netech.org>

Industrial Site: <http://www.PressureMeter.com>